A Calculation about preparation of stock solution

Stock solution for standard

Stock solution 10,000 ppm standard As(III):

: NaAsO₂(MW) = 129.91017, As (MW) = 74.9216
10,000 ppm = 10,000 mg/l or 1 g/100 ml
As = 74.9216 g, NaAsO₂ = 129.91017
As = 1 g, NaAsO₂ =
$$(1x129.91017)/74.9216$$

= 1.7340 g

So weighing 1.7340 g of NaAsO₂ and dissoving in 100 ml of deionize water, will give stock solution 10,000 ppm standard As(III).

Pipet 10 µl stock solution 10,000 ppm standard As(III)

: Standard As(III) 10,000 ppm : 10^3 ml, have As 10^4 mg Pipet standard As(III) 10 μ l (0.01 ml), have As = [(10^4 mg)x(0.01 ml)]/(10^3 ml) = 0.1 mg

So pipet standard As(III)10 μ l from stock solution 10,000 ppm standard As(III) will have As 0.1 mg.

B Sources of samples

Sample number	Sources of samples
1	Well water from Moo 2
2	Creek water from Moo 2
3	Marsh water from Moo 2
4	Creek water from Moo 2
5	Well water from Moo 12
6	Ground water from Moo 12
7	Marsh water from Moo 2
8	Well water from Moo 7

The water samples were from Ronpiboon district, Nakhon Sri Thammarat province, Thailand.

VITA

Name -

Miss Ratiya Singhon

Date of birth

June 30, 1976

Academic status

High school certificate holder from Benjamamaharaj, Ubon

Ratchathani, 1995

B.S. (Chemistry) degree holder from Naresuan University,

1999

M.S. (Chemistry) degree holder from Chiang Mai University,

2002

Paper Presented

R. Singhon, S. Liawruangrath, Th. Chuesaard and B.

Liawruangrath, Determination of aluminium(III) in

preparation by atomic absorption spectrophotometry, 26th,

Congress on Science and Technology of Thailand, 2000.