



APPENDIX

ลิขสิทธิ์มหาวิทยาลัยเชียงใหม่

Copyright© by Chiang Mai University

All rights reserved

Reagent preparation for phytochemicals screening

1. Dragendorff's reagent

Preparation of bismuth nitrate solution, 8.0 g bismuth nitrate in 12.0 mL 30% nitric acid. Dissolve 27.2 g potassium iodide in 50 mL water and put in to the bismuth nitrate solution then diluted to 100 mL with water.

2. Kraut's spray reagent

Dissolve 8.0 g bismuth nitrate in 20 mL nitric acid and 27.2 g potassium iodide in water 50 mL. mix the two solutions and allow precipitation, then take the clear solution and adjust to 100 mL with water.

3. Mayer's reagent

Two solutions; 1.36 g mercuric chloride in 60 mL water and 5.0 g potassium iodide in 10 mL water are mixed and make the volume diluted to 100 mL with water.

4. Valser's reagent

Gently dissolve 2.0 g mercuric iodide in potassium iodide solution (10.0 g KI in 80 mL water) and adjust volume to 100 mL with water.

5. Wagner's reagent

Dissolve 2.0 g potassium iodide in water, add 1.27 g crystals of iodine then shake and adjust volume to 100 mL with water.

Preparation of spray reagents

1. Anisaldehyde in sulfuric acid

Add 1 mL conc. Sulfuric acid to the mixture of 0.5 mL anisaldehyde in 50 mL glacial acetic acid. Freshly prepare before using.

2. Dragendorff's spray reagent

Solution A: add 10 mL water into the solution of 0.6 g bismuth subnitrate in 2 mL conc. hydrochloric acid.

Solution B: dissolve 6.0 g potassium iodide in 10 mL water and shake.

Mix the solution A and B then add 7 mL conc. Hydrochloric acid and 15 mL water. Adjust the volume to 400 mL with water.

C:\Documents and Settings\Administrator\My Documents\NMR-data\M2\chata\ON\2\fi-0\SKDammrong-ONC-G6.als
G6

DFILE C:\Documents and Settings\Adm
COMNT G6
DATIM 29-MAY-2009 16:33:20
OBNUC 1H
EXMOD single_pulse.exp
OBFREQ 600.00 MHz
OBSET 175.00 KHz
OBFIN 905.7 Hz
POINT 32768
FREQU 9009.0 Hz
SCANS 16
ACQTM 3.637 sec
PD 4.000 sec
PWL 3.1 us
IRNUC 1H
CTEMP 17.8 C
SLVNT CDCL3
EXREF 0.00 ppm
BF 0.01 Hz
RGAIN 17




Figure 19 The 600 MHz ^1H -NMR spectrum of compound G6P

C:\Documents and Settings\Administrator\My Documents\NMR-data\M2\chana\ON\2&f-0\SKDamrong-CNC-G6-BCM.als
G6_BCM

DFILE C:\Documents and Settings\Adm
COMNT G6 BCM
DATIM 29-MAY-2009 18:22:42
OBNUC 13C
EXMOD single_pulse_dec
OBFRQ 150.00 MHz
OBSEF 930.00 KHz
OBFIN 32.7 Hz
POINT 65536
FREQU 37735.8 Hz
SCANS 2001
ACQTM 1.737 sec
PD 1.500 sec
PW1 9.2 us
IRNUC 1H
CTEMP 18.8 C
SIVNT CDCL3
EXREF 77.00 ppm
BF 2.00 Hz
RGAIN 27




Figure 20 The 150 MHz ^{13}C -NMR spectrum of compound G6P




Figure 21 The 600 MHz ^1H -NMR spectrum of compound P6bP




Figure 22 The 600 MHz $^1\text{H-NMR}$ spectrum of compound P6fP

C:\Documents and Settings\Administrator\My Documents\NMR-data\M2\ohara\ON\2ø-0\SKDammrong-ONC-P6fP-BCM.als
P6fP BCM

C:\Documents and Settings\Adm
COMNT P6fP BCM
DATIM 29-MAY-2009 16:13:35
OBNUC 13C
EXMOD single pulse dec
OBFRQ 150.00 MHz
OBSET 930.00 KHz
OBFIN 32.7 Hz
POINT 65536
FREQU 37735.8 Hz
SCANS 889
ACQTM 1.737 sec
PD 1.500 sec
PW1 9.2 us
IRNUC 1H
CTEMP 19.1 c
SLVNT CDCL3
EXREF 77.00 ppm
BF 2.00 Hz
RGAIN 27




Figure 23 The 150 MHz ^{13}C -NMR spectrum of compound P6fP




Figure 24 The 600 MHz ¹H-NMR spectrum of compound H4d1P



```

----- ACQUISITION PARAMETERS -----
File Name      = SDA000009-ONE-H4D1P
Author        = JEOL LTD.
Sample ID     = SDA000009-ONE-H4D1P
Content       = ONE-H4D1-PCM
Creation Date  = 7-APR-2009 17:25:28
Revision Date = 7-APR-2009 17:30:21
Spec Site     = ECP600

Spec Type      = DELTA_NMR
Data Format     = 1D_COMPLEX
Dimensions     = X
Dir Title      = 13C
Dir Name       = 65536
Acq Units      = Hz
Acq Unit Start Time = 7-APR-2009 17:25:29
Delay of Start = 1[s]
Digital Filter = FALSF
End Time       = 7-APR-2009 17:34:07
Field Strength = 14.05636928[T]
Filter Mode    = BUTTERWORTH
Filter Width   = 18.8641788[MHz]
Irr_Solvent    = CDCl3
Irr_Freq       = 600.17530548[MHz]
Irr_Noise      = NAUGH
Irr_Offset     = 8.0[ppm]
Irr_Width      = 58[us]
Iterations     = 1
Local Time     = 7-APR-2009 17:25:28
Obs_Noise      = NAUGH
Proc_Width     = 3188
Recvr_Gain     = 29
Relaxation_Delay = 1.5[s]
Scans          = 5841
Solvent        = CHLOROFORM-D
Spin_Set       = 14[Hz]
Spin_Lock_90   = 25[us]
Spin_Lock_Attn = 23[dB]
x90_Awt        = 18.5[us]
X_Acq_Duration = 1.736704[s]
X_Angle        = 45[deg]
X_Domain       = 13C
X_Offset       = 150.9300327[MHz]
X_Preamplifier = 65536
X_Preamplifier Gain = 4
X_Pulse        = 9.25[us]
X_Resolution   = 0.57581215[Hz]
X_Sweep        = 37.72584906[MHz]
T1_Rho        = 10[us]
T1_Noise      = NAUGH
T2_Rho        = 1[us]
T2_Noise      = NAUGH
Coa_Width     = 1[us]
  
```




Figure 25 The 150 MHz ¹³C-NMR spectrum of compound H4d1P




Figure 26 The 600 MHz ^1H -NMR spectrum of compound I8eP




Figure 26 The 150 MHz ¹³C-NMR spectrum of compound I8eP


VITA

NAME Mrs. NICHAKAN PHOOWIANG

DATE OF BIRTH May 22, 1977

PLACE OF BIRTH CHIANG MAI, THAILAND

INSTITUTION ATTENDED BACHELOR OF EDUCATION OF SCIENCE
RAJABHAT CHIANG MAI INSTITUTE



ลิขสิทธิ์มหาวิทยาลัยเชียงใหม่

Copyright© by Chiang Mai University

All rights reserved